

DISINFECTANTS

VIRUSOLVE® + WIPES

DISINFECTION & CLEANING FOR ANY SURFACE & NON INVASIVE MEDICAL DEVICES

Virusolve®+ does not contain any hazardous aldehydes or chlorine generating components. Virusolve®+ is highly effective as a cleaner effective against mycobacteria, bacteria, fungi and viruses. The product is designed to eliminate the risk of cross-infection from TB, viruses, fungi and bacteria.

Virusolve®+ does not simply render the microorganisms inert but has been demonstrated to kill these by disrupting the RNA of the microorganisms.

It is safe to use on any surface and is odourless and colourless.



Primary features attributable to Virusolve®+ are:

- Disinfection and cleaning.
- Does not contain chlorine or aldehyde compounds.
- Effective over a wide range of pathogens including mycobacteria, viruses, bacteria, fungi, yeasts and moulds, including MRSA, Hepatitis B and HIV.
- Effective against the difficult non-enveloped type viruses such as Polio, Adenovirus and Norovirus.
- Non-fixative.
- Non-corrosive - will not damage delicate coatings.
- No activation required.
- Cost effective.
- No fumes.

Packaging

- Ready to Use
- 200 soft pack
- 225 tub

DISINFECTANTS

Application

Virusolve®+ wipes are ideal for removing biohazards in hospitals, health centres, dentists, vets, nursing homes, laboratories, food preparation areas and any environment where the risk of cross-infection needs to be eliminated, e.g. cruise ships, disinfecting hairdressing equipment etc.

Virusolve®+ wipes are ideal for the cleaning and disinfection of hard surfaces, walls, floors and glasswork.

Directions for Use

Pre-impregnated Virusolve®+ wipes are supplied in containers and are ready to use direct from the container.

General Surface Cleaning and Healthcare Environment:

Take a wipe and wipe gently over the whole area to be cleaned. Allow a minimum contact time of 1 minute before wiping clean to remove grease, blood and organic debris. Removal of stubborn deposits may be assisted by the use of a stiff bristle brush. Rinse treated surface with clean water if required.

Note: Food preparation surfaces should be allowed to dry before use.

For Disinfection of medical, surgical and dental instruments and equipment (Manual procedure):

Wipe all surfaces of the instrument (preferably disassembled, where applicable) using a pre-impregnated Virusolve®+ Wipe and allow the fluid to remain in contact for at least 5 minutes.

Examine the instrument and wipe or brush clean to remove debris that remains.

Rinse treated instrument in sterile, cold water and allow to dry.

Replace used Virusolve®+ wipes frequently or when gross soiling of the wipe becomes evident.

VIRUSOLVE®+ WIPES DISINFECTION & CLEANING FOR ANY SURFACE & NON INVASIVE MEDICAL DEVICES

Virusolve®+ Properties

| Active Against | Standards | Contact Time |
|----------------|--|-------------------------------------|
| Mycobacteria | EN14563 Mycobacterium terrae Mycobacterium avium | 15 Minutes 15 Minutes |
| Bacteria | EN1276, EN13697, EN13727, EN14561, EN 16615 | 1 Minute 5 Minutes |
| Viruses | EN 14476 H1N1 & Noro | 5 Minutes |
| Fungi / Yeasts | EN1275, EN1650, EN13697, EN16615 EN14562 | 1 Minute 5 Minutes 15 Minutes |

Active ingredients: Alkyl dimethyl benzyl ammonium chloride, Polymeric biguanide hydrochloride.



Amity Limited, Libra House, West Street,
Worsbrough Dale, Barnsley, South Yorkshire,
S70 5PG, UK

T: +44 (0) 1226 770787

F: +44 (0) 1226 770757

E: sales@amityinternational.com

WWW.AIHC.CO.UK

© Amity Limited all rights reserved 2018

DISTRIBUTED BY: