

## DISINFECTANTS

# VIRUSOLVE® + TRIGGER SPRAY

## HIGH LEVEL DISINFECTION & CLEANING FOR MEDICAL DEVICES AND ANY SURFACE

Virusolve®+ is a patented proven HLD (High Level Disinfectant). The formulation does not contain any hazardous aldehydes or chlorine generating components.

Virusolve®+ is highly effective as a cleaner and terminal HLD against bacterial spores, mycobacteria, bacteria, fungi and viruses. The product is designed to eliminate the risk of cross-infection from spores, TB, viruses, fungi and bacteria.

Virusolve®+ does not simply render the microorganisms inert but has been demonstrated to kill these by disrupting the RNA of the microorganisms. It is safe to use on any surface and is odourless and colourless.

### Primary features attributable to Virusolve®+ are:

- High Level Disinfection and cleaning.
- Does not contain chlorine or aldehyde compounds.
- Effective against bacterial spores e.g. Clostridium difficile (C dif), Bacillus cereus etc.
- Effective over a wide range of pathogens including mycobacteria, viruses, bacteria, fungi, yeasts and moulds, including MRSA, Hepatitis B and HIV.
- Effective against the difficult non-enveloped type viruses such as Polio, Adenovirus and Norovirus.
- Non-fixative.
- Non-corrosive - will not damage delicate coatings.
- No activation required.
- Cost effective.
- No fumes.
- Tested to EN14885 for disinfection properties



### Packaging

#### Ready to use

- 750ml trigger spray
- 500ml trigger spray



**Guaranteed protection and performance,  
giving total infection control - safely**

## DISINFECTANTS

### Application

General Surface Cleaning and Healthcare Environment:

For High Level Disinfection: Use Virusolve®+ trigger spray

Spray Virusolve®+ solution onto a cloth/wipe and apply to the surface being treated. Wipe over the whole area.

Allow a minimum contact time of 1 minute before wiping clean or brushing to remove grease, blood and organic debris. Removal of stubborn deposits may be assisted by the use of a stiff bristle brush. Rinse treated surface with clean water if required.

Virusolve®+ solution has demonstrated efficacy in the presence of organic soil contamination and microbiological burden during reuse. The active constituents' concentration of Virusolve®+ solution during its use-life must be verified by the Virusolve®+ Concentration Test Kit or Test Strips.



## VIRUSOLVE®+ TRIGGER SPRAY

HIGH LEVEL DISINFECTION & CLEANING FOR  
MEDICAL DEVICES AND ANY SURFACE

### Virusolve®+ Properties

Active Against	Standards	Contact Time
Spores	EN13704 EN14347	2 Minutes 1 Minute
Mycobacteria	EN14348 Mycobacterium terrae Mycobacterium tuberculosis EN14563 Mycobacterium terrae Mycobacterium avium	1 Minute 1 Minute 5 Minutes 5 Minutes
Bacteria	EN1276, EN13727, EN13697, EN14561	1 Minute 5 Minutes
Viruses	EN 14476 Polio, Adeno & Noro HIV, Hep B & C	5 Minutes 1 Minute
Fungi / Yeasts	EN1275, EN1650 EN13624, EN14562	1 Minute 5 Minutes

Active ingredients: Bis (3-aminopropyl) dodecylamine.  
Didecyldimethyl ammonium Chloride



Amity Limited, Libra House, West Street,  
Worsbrough Dale, Barnsley, South Yorkshire,  
S70 5PG, UK

T: +44 (0) 1226 770787

F: +44 (0) 1226 770757

E: sales@amityinternational.com

[WWW.AIHC.CO.UK](http://WWW.AIHC.CO.UK)

© Amity Limited all rights reserved 2018

DISTRIBUTED BY: