

DISINFECTANTS

AMITY PAA 15


SINGLE PROCESS PERACETIC ACID, NO ACTIVATOR REQUIRED

Amity PAA 15 is a powerful peracetic acid based high level disinfectant / sterilant concentrate for use in Automated Endoscope Reprocessors (AERs). It is an acidic disinfection concentrate for flexible, rigid endoscopes and medical devices.

Amity PAA 15 is registered as Class IIB Medical Device

Amity PAA 15 is a HLD fluid effective against, spores, mycobacteria, viruses, bacteria and fungi when used at a concentration of 1%.

Devices and endoscope parts should be cleaned prior to disinfection with Amity PAA 15 using suitable cleaners such as Viruzyme family or Virudet N.



Primary features attributable to Amity PAA 15 are:

- Powerful peracetic acid based high level disinfectant / sterilant (sporicidal, mycobactericidal, virucidal, bactericidal and fungicidal).
- Recommended dilution 10ml/L or 1% in RO water.
- Single component & does not require activator compounds.
- Contains anti-corrosion additives.
- Rapid & wide ranging activity including spores.
- HLD claim: 5 minutes at 35°C (10 minutes at 20°C).
- Tested to EN 14885 for disinfectant properties.
- Tested to EN 15883-4.
- Contains no aldehydes and formaldehyde.

Packaging

- 5 litre containers



Amity PAA 15 - A Peracetic Acid based High Level Disinfectant (HLD) / Sterilant Fluid for Medical Devices & Instruments

DISINFECTANTS

Application

For disinfection of Medical Instruments and Endoscopic Equipment

Stage 1:

Wash medical device carefully in an enzymatic cleaner (Viruzyme Family or Virudet N) and rinse.

Stage 2:

Amity PAA 15 should be dosed at 10ml/L into AER.

Phase	Time	Temp.	Solution	Dosage
Wash	5 min	25-55 °C	Viruzyme	5ml/L (0.5%)
Rinse	3 min	-	-	-
Disinfection	5 min 10 min	35 °C 20 °C	Amity PAA15	10ml/L (1.0%)
Rinse	3 min	25 °C	-	-
Rinse	3 min	25 °C	-	-
Dry	-	-	-	-

Drying time and temp according to the machine manufacturer

AMITY PAA15

SINGLE PROCESS PERACETIC ACID, NO ACTIVATOR REQUIRED

Amity PAA 15 Properties

Active Against	Standards	Contact Time
Spores	EN17126	5 Minutes
Mycobacteria	EN14348	5 Minutes
	Mycobacterium terrae	5 Minutes
	Mycobacterium avium	5 Minutes
	EN14563	5 Minutes
Mycobacterium terrae		5 Minutes
	Mycobacterium avium	5 Minutes
Bacteria	EN13727,	5 Minutes
	EN14561	5 Minutes
Viruses	EN 14476	5 Minutes
	Polio & Adeno	5 Minutes
	HIV, Hep B & C	
Fungi / Yeasts	EN13624	5 Minutes
	EN14562	5 Minutes

Active ingredients: Peracetic Acid, Hydrogen peroxide.



UK
CA
0120

CE
1639

Amity Limited, Libra House, West Street,
Worsbrough Dale, Barnsley, South Yorkshire,
S70 5PG, UK

T: +44 (0) 1226 770787

F: +44 (0) 1226 770757

E: sales@amityinternational.com

WWW.AIHC.CO.UK

© Amity Limited all rights reserved 2023

DISTRIBUTED BY: