



Viruzyme[®] Ve

A Protolytic Enzymatic Biological Cleaner for Manual & Automated Cleaning

Technical Information

Description & Application:

Viruzyme Ve is a protolytic enzymatic detergent for use in washer/disinfectors, endoscope Reprocessors & decontamination systems and manual decontamination.

It is a slightly alkaline enzymatic detergent formulated to wash instruments that may be damaged by highly alkaline chemistries. It provides superior cleaning action without foaming or causing corrosion on instruments or washer/disinfectors and has a wide material compatibility. It is compatible with stainless steel, aluminium, anodised aluminium, chromium-plated material, synthetic materials, rubber, latex, glass and ceramics.

Viruzyme Ve is particularly effective for the removal of protein and organic contaminants from sensitive surfaces, including dried on blood and tissue.

Viruzyme Ve is registered as a Class I medical Device

It is suitable to use on any surface and has an effective cleaning action during manual, ultrasonic, automated washer/disinfector and AER equipment.

It is normally used in the wash phase after the pre-rinse stage of automated washer/disinfector equipment, and does not interfere with the cleaning mechanism in ultrasonic and washer/disinfector equipment.

Primary features attributable to Viruzyme Ve are:

- Biological cleaner suitable for removal of proteinaceous and organic contaminants from sensitive instrument materials;
- Low foaming;
- Suitable for manual, ultrasonic, washer/disinfector & AER equipment
- Uses a blend of five enzymes, protease, amylase, lipase, cellulase and mannanase and synergistic surfactants that have added safety:
 - Low toxicity
 - Pleasant odour
 - Reduced instrument damage and repair
 - Excellent material compatibility
 - Non-fixative
 - Non-corrosive, will not damage delicate coatings



Certificate No. GB06/69741



Certificate No. GB06/69740



**Infection
Controlled**

- Highly effective cleaner action

Applications

Viruzyme Ve is ideal for removing bioburden from medical instruments, endoscopes and other equipment including those with delicate or sensitive surfaces that are used in hospitals, health centres, dentists and laboratories.

Compatibility

Viruzyme Ve is compatible with a wide range of fabrics, rubbers, plastics and hard surface materials and will not attack metal surfaces when used in accordance with standard instructions for use.

Health and Safety

Viruzyme Ve is biodegradable.

Viruzyme Ve is certified as being non-hazardous when diluted and used according to the specified user directions.

The organisms used in the formulation are naturally occurring and are known to be non-pathogenic to humans and animals.

See Safety Data Sheets for details.

Typical Properties

Appearance	Pale blue liquid
Odour	Low odour
pH,	10.5 units (Neat solution)
Density, g/cm ³ @ 20°C	1.07 (neat)
Solubility in water	Complete
Flash Point (Abel closed cup)	None

Shelf Life: Concentrate – 2 years in unopened original containers when stored between 5°C and 25°C out of direct sunlight.

Directions for Use

Viruzyme Ve is supplied as a concentrated product that requires dilution prior to use.

Application:

For Manual Cleaning of Instruments and Equipment:

Use Viruzyme Ve product at a dilution of 6 to 10 ml/L in water (i.e. 0.6 to 1.0% v/v)

Greater strengths up to 40ml/L (4% v/v) may be used where high levels of contaminants or where difficult or hazardous contaminants may be present.

Soak parts for between 5 and 10 minutes at 20°C to 30°C in the diluted Viruzyme Ve solution.

(Note: allow longer periods – up to 20 minutes when dealing with dried on proteinaceous deposits or if low solution temperatures are present).

Finally, rinse with clean ware and transfer to next cleaning stage.

For Ultrasonic Cleaning of Instruments and Equipment:

Use Viruzyme Ve product at a dilution of 6 to 10 ml/L in water (i.e. 0.6 to 1.0% v/v)

Greater strengths up to 40ml/L (4% v/v) may be used where high levels of contaminants or where difficult or hazardous contaminants may be present.

Soak parts for between 5 and 10 minutes at 20°C to 50°C in the diluted Viruzyme Ve solution. (Note: allow longer periods – up to 20 minutes when dealing with dried on proteinaceous deposits or if low solution temperatures are present).

Finally, brush surfaces to remove any gross debris then rinse thoroughly with clean water and transfer to next cleaning stage (e.g. washer/disinfector, ultrasonic or AER equipment).

For Automated Cleaning of Instruments and Equipment:

Injectors should be adjusted to deliver between 6 to 10 ml/L in water (i.e. 0.6 to 1.0% v/v) dependant on the expected soil loading of the contaminated equipment and the temperature of the water.

A recommended automatic cleaning program is as follows:

Phase	Water type	Time	Temp	Detergent	Dosage
Pre-wash	Cold tap water	3 min	-	-	-
Wash	Hot tap water	5 to10 min	≤55°C	Viruzyme Ve	6 - 10ml / L (0.6 – 1.0%)
Rinse	Hot water	2 min	-	-	-
Rinse	DI or RO water	2 min	-	-	-

Viruzyme Ve is not a disinfectant and parts should proceed to disinfection or sterilisation by appropriate methods immediately following cleaning.

A recommended AER cleaning program is as follows:

Phase	Water type	Time	Temp	Detergent	Dosage
Pre-wash	Cold tap water	3 mins	-	-	-
Wash	Hot tap water	5 to10 mins	≤55°C	Viruzyme Ve	6 - 10ml / L (0.6 – 1.0%)
Rinse	Water	2 mins	-	-	-
Disinfection (HLD)	DI or RO water	5 to 10 mins	20°C 35°C 35°C 35°C 25°C	Virusolve+ EDS Auto Amity PAA 15 Amity PAA 5 System Amity PAA 10 Amity PAA RTU Amity OPA	25ml/L 10ml/L 20ml/L part A/B 15ml/L Neat (no dilution) Neat (no dilution)
Rinse	DI or RO water	2 mins	-	-	-
Rinse	DI or RO water	2 mins	-	-	-

General:

Initial rinse temperature should be below 35°C otherwise the temperature may bake on protein contaminants rendering them more difficult to remove.

Enzyme activity increases with increased water temperature up to a maximum temperature of 55°C.

A policy by the practice to carefully monitor the above will help to reduce any potential problems.

Do not add other cleaning agents to the Viruzyme Ve solution as they can seriously impair performance of the Viruzyme Ve cleaner.

Disposal of spent solutions:

Flush to drain with copious water or soak up onto inert material and dispose of with clinical waste
Product is biodegradable under OECD conditions.

CONTACT DETAILS:

For UK and Rest of the World:

Amity International,
Libra House, West Street
Worsbrough Dale,
BARNSELY
S YORKS, S70 5PG,
ENGLAND

Tel: +44 (0) 1226 247660
Fax: +44 (0) 1226 291936

E-Mail: sales@amityinternational.com

Website: www.amityinternational.com

For North America:

Amity International,
PO Box 5254,
1704 Denver Road,
ANDERSON,
SOUTH CAROLINA,
SC29623, USA.

Tel: 864 622 2233
Fax: 864 622 2234

For any further information, please contact your distributor or Amity.

In the event of any technical queries, please contact:
Mr. Ram Singh at the UK/ROW address or by e-mail to:

rsingh@amityinternational.com

Packaging Details:

Viruzyme Ve is available in the following packaging.

- 5 L container
- 10 L container (available on request)
- 25 L container (available on request)
- 200 L drum (available on request)